1A400

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1A400

Assembling Instructions


SAFETY INSTRUCTIONS FOR THE PATIENT

- The cosmetic shell is guaranteed for 6 months. It must be replaced if it gets torn or if there is some slack between the structure and the cosmetic shell.
- Do not let foreign bodies or dust inside the foot because they could lead to noises or early abrasions. They must be removed from the cosmetic shell.
- If the foot has been dipped or sprinkled, empty the remaining liquid out and dry the foot rapidly. In case of metal component corrosion, have your prosthesis controlled.
- Use temperature : -20° / $+60^{\circ}$ C
- An inappropriate use of the foot (top level sport, frequent overload, repeated shocks) can damage components and be dangerous for you. If you discover an abnormal behaviour of your foot or changes in your foot features, contact your prosthetist.

1- Description and functions

Dyna C foot includes three main components. Two composite parts : a C-form malleolus that allows energy storage and return, and a sole providing flexibility on the sides, which ensures adaptation to the ground. A polyurethane cosmetic shell which transfers the mechanical stresses between the ground or the shoe and the composite structure. Thanks to these components, the foot behaviour is flexible and dynamic. For patients who are active, or very active (children for example), or reach the weight and activity limits of the module selection guide, a stop can be used to limit warping of the C-malleolus. This stop comes with the foot and must be installed by the prosthetist.

2- Specification

Dyna C foot is designed exclusively for lower limb prostheses. It is appropriate for active patients up to 125 Kg, or very active patients up to 105 kg.


Dyna C foot is designed for a 10 mm heel. It includes a male pyramid that allows for the use of standard connectors available on the market. Thanks to these connectors, the heel height can be easily modified from 0 to 20 mm.

Module selection guide


Activity	Weight				
	50 to 60	60 to 75	75 to 90	90 to 105	105 to 125
Moderate	1	1	2	3	4
Normal	1	2	3	4	5
High	2	3	4	5	N/A

3- Assembling (without the patient)


A) Gluing of the stop (if necessary): Remove the foot from its cosmetic shell, then take the sock off. Degrease the faces that will be glued, let them dry, then apply epoxy glue on both faces.

Put the stop in contact with the nut (1), then position it by pushing it upwards (2) and towards the nut (3). Let the glue dry, then put the sock and the cosmetic shell again.

B) Assembling: Place under the foot a wedge of the same height as the shoe heel. Fix the female pyramid connector to the foot taking care that its axis be perpendicular to the ground (fig1).

4- Alignment (with the patient): Place the foot so that the load line goes through the posterior third of the foot (fig. 2)



Fig. 2

Notes


- Increasing equinus provides a more dynamic alignment. It increases energy return at step completion.
- If the foot seems too long during transfer to the forefoot or if heel support is too flexible, you have to position the foot behind the socket.
- If heel support seems too hard, or if the foot seems too short during transfer to the forefoot, you have to position the foot ahead of the socket.

5- Finishing: To avoid rubbing noises, pull the sock up to the ankle connector and fix it with a fine string for example. The malleolar clips that ensure connection to the cosmetic foam allow to dismantle the foot easily and rapidly.

	Foot size (cm)	22 to 25	26 to 27	28 to 29
	Product Nr	1A40013	1A40014	1A40015

Installation: Cut the foam to the desired length, position the foam on the foot and draw the malleolus contour. Inside the foam mill a space that will receive the selected connection part and the foot part protruding from the shell, and at the same time take care that thickness be sufficient to allow for gluing. Position the malleolar clip carefully and clip it into the cosmesis. Clean contact surfaces with isopropyl alcohol and glue the malleolar clip on the foam with XC090 neoprene adhesive. After drying, finish the outer shape as usually.

6- Maintenance: The cosmetic shell, the protective sock and the malleolar clip are the only foot components that may be worn out and may need to be replaced. The cosmesis can be removed and put again with 1A40019 insert horn.

Spare cosmesis 1A40020-D22 including the cosmetic shell, the sock and the malleolar clip.
Side* Size in cm

7- Recycling : The different components of the foot are special wastes and must be treated accordingly.

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